

Baby Boutique Online SA - Showroom : Glen Marais,
Kempton Park, JHB (by appointment only)
Factory : Cilvale AH, Bapsfontein, JHB (by appointment
only)

Phone: 0761588448 Rene - Email:
enquiries@babyboutiqueonline.co.za

Lumilia Lux 2in1 IPL Hair Removal & Skin Rejuvenation Handset - Pink Rose

Price: ~~R3 299.00~~ R2 399.00

Description

IPL is suitable for all but the darkest skin tones. This is because the light energy may be more attracted to the pigments in the skin rather than the melanin in the hair. Refer to the skin tone chart on below to see whether or not IPL is suitable for you.

The Lumilia Lux™ works for all but the darkest skin tones. Kindly refer to the Skin & Hair tone guide below to see if your skin type is a fit.



The Lumilia Lux™ is not effective if your hair colour is light blonde, grey, red or white. Check the hair colour chart below to guide you.



Out with the old and in with the new. Meet the all-new **Lumilia Lux 2-in-1**, boasting with advanced **@Home IPL Laser Hair Removal Technology & a new Skin Rejuvenation mode**. Why settle for less if you can **invest in the best**.

The **first and only** device to offer **8 different energy levels** with an **easy to access color LCD screen** to switch between **Hair Removal and Skin Rejuvenation modes**. An all-new **luxury design, more powerful** than any of its predecessors, including a **quieter** and more efficient **dual-fan** to accommodate the **larger** light window. All of this to give you more **FLASH** for your investment!

Designed to deliver **long-lasting visible results in just 2-3 weeks**, full-body

treatments are done in less than 30 minutes, as little as once per week.

Exceptional, Innovative, and Effective.

Works on any body part including the **face** and **Brazilian area**. Your handset will utilize clinical-grade technology that is **certified** for **safe and effective** at-home hair removal.

Get **long-lasting results** at a **fraction of the cost of clinical laser treatment**.

There are **no hidden costs** and **no refills or replacements** needed, ever! Invest in a new Sleek and Silky smooth you!

?? HAIR REMOVAL & SKIN REJUVENATION - The first and only device to offer a Hair Removal and Skin Rejuvenation mode. Skin Rejuvenation can assist to reduce and improve the appearance of hyperpigmentation, wrinkles, reduce skin redness, acne scars, sunspots, boost collagen, etc.

?? NO PAIN - IPL Laser feels like a warm sensation/zap on the skin. It is nothing compared to Razor Bumps and Burns and not even to mention those painful visits to the beauty clinic!

?? SAFE - It is clinically proven to be safe for your skin. You can even use it on sensitive parts like your face, armpits, and bikini area.

?? 1-8 INTENSITY LEVELS -Choose between 8 different energy levels. Select the intensity level that suits you best.

?? 20+ Years 800,000 Flashes - Lumilia Lux IPL Handset has 800,000 Light Pulses which will last you a lifetime

FAQ's

What is IPL Laser Hair Removal and how does it work

IPL (Intense Pulsed Light) is a non-invasive form of light therapy used for various dermatological procedures including long-term/ permanent hair removal. The device emits a broad spectrum of light which is absorbed by the pigment melanin located in the hair follicle. When the pigment in the hair follicle absorbs the IPL light, the light energy is converted into heat, resulting in thermal damage to the hair follicle. This results in delayed regrowth of the hair or permanent destruction of the hair follicle. Treated hairs naturally fall out over the course of a few weeks.

How strong is this product and what is the energy output?

Our device is built with advanced features that other devices won't have: • LED Display Screen to visibly view the device settings • Latest clinical grade technology • Skin Rejuvenation Setting, which improves collagen production, skin tone & texture, acne scars, skin redness, spider veins, hyperpigmentation, etc. • 3.5cm large treatment window • Maximum energy output of 15j for stronger and more effective treatments • Wavelength spectrum of 550-1200nm which ensures entire hair follicle is treated • Built in Integrated Safety Sensor to prevent unintentional flashing • Glide and Auto Flash Mode for treating large areas fast • 8 Adjustable Light Settings to suite different skin tones

Which areas of my body can I use this on?

You can use your handset on your entire body, including your Brazilian and face (just be sure to not get too close to your eyes).

Do I need to shave before using this product?

Yes, you should shave the desired areas before using your handset. A close to skin shave is best recommended. The goal is for the light energy to reach the hair follicle as quickly as possible, so shaving beforehand is essential to achieve the best results. The hair follicle does need to be in tact for IPL to be effective. This means no waxing, plucking, threading, epilating etc.

Which hair removal methods can I use between treatments?

Between treatments you should only shave when required. You should never wax, pluck, or epilate as these hair removal methods removes the entire root, which is what absorbs the light during the treatment.

Are there any reoccurring costs?

No. This product is a once-off purchase. There is absolutely no need for replacements or refills.

Does IPL work for all skin tones?

IPL is suitable for all but the darkest skin tones. This is because the light energy may be more attracted to the pigments in the skin rather than the melanin in the hair. Refer to the skin tone chart on the How To Use page to see whether or not IPL is suitable for you.

How long will this device last?

Our upgraded handsets have over 800,000 flashes, which will last well over 20 years if used as directed. The device comes with a 1-year warranty and a 90-day money-back guarantee.

What is the Skin Rejuvenation setting?

This device also includes a unique Skin Rejuvenation Setting at the Level 1 Light Intensity which can be used weekly. The device releases a pulse of filtered light that stimulate the fibroblast cells within the skin that produce collagen and elastin. By stimulating collagen production, this results in an improved and more youthful skin tone and texture. It is effective in controlling fine lines, acne and sun damage.

Can I use my Lumilia IPL Handset before or after tanning?

As IPL works with your natural skin tone, adding artificial tan creates a layer over your skin which doesn't allow the device to penetrate as effectively. For a successful treatment please wait 2 weeks after sun tanning to use your Handset and please also wait at least 48 hours to tan after using your handset. Artificial tan should also not be present on your skin whilst using your handset but can be applied after.

How long is the warranty for the Lumilia IPL Handset?

The Lumilia Handset comes with a 12-month warranty so if your device stops flashing,

we will replace the device for you free of charge. Please visit our Shipping and returns page for terms and conditions.

How often are treatments?

Use your Lumilia Handset every 2-3 days, for your first 4-5 treatments. Once the 4 to 5 treatments are completed, we recommend touch ups every 2nd week. IPL is not a 100% permanent, so you may reach the stage where very minimal hairs grow back. Skin Rejuvenation can be used 3 times per week.